



IS-III active
S-narrow




Super Narrow for Anterior Mandible (Narrow Ridge)



Diameter	Length	Product name
Ø3.2	8.5mm	IS3SN3008AP
	10.0mm	IS3SN3010AP
	11.5mm	IS3SN3011AP
	13.0mm	IS3SN3013AP

* Cover screw not included.

IS-III active S-Narrow Features




S.L.A. Surface

The new S.L.A Surface with 40% greater surface area and 50% more cell adhesion promotes osseointegration



Straight Body

Easy to adjust implantation depth



Deep & Wide Pitch

Optimum Pitch for osseointegration (0.9Pitch)




Wide Cutting Edge

Wide cutting edge and enlarged surface area enhances initial fixation and offers clinicians more stable implant placement




• M2.0 Connection

- One Piece Type
- Cover Screw
- Healing Abutment



• Internal Conical Seal 11° & 2/1 Hex



• M1.6 Connection

- Two Piece Type
- PickCap Impression Coping
- Temporary Abutment
- Cement Abutment
- Angled Abutment



S-narrow Prosthesis

- IS-III active S-narrow
- Lab Analog
- Cover Screw
- Healing Abutment
- PickCap Impression Coping
- Temporary Abutment
- Cement Abutment
- Angled Abutment

