


Configuration	A3	A3L	A3S
Electrolytic Oxidized Disinfection System	✓	△	△
Infrared & Gravity Dual Induction Water Supply	✓	△	△
Thermostatic water system for dentist table	✓	△	△
Air Regulator	AirTAC	△	△
Timer	✓	-	-
Three memory positions	✓	-	-
Air & Water Tube	SMC	Polyether	Polyether
Dental light	Rolight E	Rolight S	Rolight S
Chair	Breathable seamless microfiber leather	Breathable seamless microfiber leather	Breathable seamless microfiber leather
	Silicone headrest cover	△	△
	Armrest holster	-	-
Built-in scaler (Digital display)	△	-	-
Built-in micromotor (Digital display)	△	-	-
Built-in micromotor (Bien Air, Digital display)	△	-	-
Dynamic atmosphere light	✓	✓	✓
Self-balance position apparatus arm	✓	✓	✓
Handpiece Holder With 4 Adjustable Angles	✓	✓	✓
Anti-collision system for Assistant table	✓	✓	△
Medical-grade color LCD screen	✓	✓	△
Smart Power On Self Test	✓	✓	△
Error Code Display	✓	✓	△
Soft start/stop system	✓	✓	✓
One-key intelligent drainage	✓	✓	✓
Ergonomically shaped Dentist stool	✓	✓	✓
Anti-suction	✓	✓	✓
Silent Round Pedal	✓	✓	✓
Four-Way Foot Control	✓	✓	✓
TIMOTION MOTOR	✓	✓	Jiechang
Stainless Steel 3-way syringe	✓	✓	△
Stainless Steel Spiral Handpiece Tube (Short)	✓	✓	✓
Nickel plated valve	✓	✓	✓
Ceramic Filter Cover	✓	✓	✓
Matte Dust-Proof Foot Pad	✓	✓	✓
5 IN1 Multifunction Tissue box	✓	△	△

Standard“✓” Optional“△” Non-optional“-”

KLT-6220 A3 Dental Unit

Grand Launch Roson Flagship Model

Made in Foshan , Loved by dentist



FOSHAN ROSON MEDICAL INSTRUMENTS CO., LTD

Add: NO.9 Henggui Mid-Road, Lianhe, Industrial Zone, Luocun, Nanhai, Foshan, Guangdong Province.

Tel: +86-757-86407037
Fax: +86-757-86407038

Web: www.fsrosen.com
E-mail: sales001@fsrosen.com





How to
develop a
great dental unit

THINKING



Foshan City
Dental Equipment Industry Base
Keep pace with the trend and future.

Roson
Dental Industry Leader
Leading the Innovation; Building the Quality.

2005年
18 years

120+
Countries

1,000,000
Users

Company Profile

Foshan Roson Medical Instruments Co., Ltd., founded in 2005, is located in Foshan---The city of manufacturing. Covering an area of more than 10000 square meters, Roson is a professional manufacturer in R&D, production and sales for dental equipment, dental chair and air compressor are Roson's main products. With 18 years' development, Roson is CE certified and has business coverage of more than 100 countries, our goal is to provide customers with the best products and services.

A3 For The Best Dentist

Foshan origin, worldwide application!

Honor Qualification





Roson Flagship Model Grand Launch


New
Appearance

Safe
Water system

Advanced
Experience

Extraordinary
Quality

KLT-6220 **A3**
Dental Unit



New Appearance
Great Innovation



Hexagonal Design



—
New looking with hexagonal element
Advanced jet black front cover
Dynamic atmosphere light
Tech-fashion appearance

Self-developed Dental light



Yellow & White

Yellow

White

—
Smiling Face design
Philips lamp beads
Stable illuminant, No RF interference, Eye protection
Surgery and Conventional treatment application
Comfortable for each treatment

INNOVATION



—
Mismatch Design
Integration of modern aesthetics and manufacturing technology
Leading the modern fashion

Enjoyable



—
Leather upholstery
Ultra-Breathable; Cosy; Environmental-friendly
As comfortable as sitting on a throne


Worry-free water Full-time disinfection

Efficient disinfection; Safety priority

EOW---Electrolytic Oxidized Disinfection System

Innovation of precedent

Water disinfection before oral usage



FOR Dentists One-Key Switch & Automatic disinfection

FOR Patients Non-toxic & No stimulation; Convenient & Safe


FOR Dental chair Corrosion resistance; Long serving lifetime

EOW-TECH


EOW-TECH Great Innovation

- Efficient disinfection**
Disinfection: $\geq 99.9999\%$
- Biofilm Inhibition**
Prevention of cross infection
- Non-toxic & No stimulation**
Without chlorine disinfectant

- Water quality adaptation**
Electrolyte adaptation
Meet different water quality
- Environmentally-friendly**
No harmful residues
- Anti-corrosion to pipeline**
No chloridion
No corrosion




Six test reports



Excellent performance
New experience



Handpiece Holder
with 4 adjustable angles (30-80°)



Convenience to pick and place, more efficient.

Thermostatic **water system**



Warm water for 3-way syringe & handpiece
One-key switch, Winter away

Anti-collision system
for Assistant table



Automatically pause and return to a section when is blocked
Easy to remove obstacles
Enable safe operations

One-key intelligent drainage




One-key press, chair up to the highest position
5 mins automatically flushed of spittoon
No extra operation, ease you on daily cleaning

Ergonomically shaped
Dentist stool

Eight-way control, easy to master
Ergonomic chair, take care dentist.

Anti-suction



Prevent water backflow
For Patient' s safety!

Medical-grade color LCD screen



Digital display/HD interface/Intelligent Timing
Built-in micromotor (optional); Built-in scaler (optional)
Synchronously display the working status and error code
Intelligent Timing for meticulous treatment

Infrared & Gravity Dual Induction Water Supply



Hidden handle
Smart dual induction system
Comfortable and time saving

Patients experience
Automatically sense and accurately measure the water every time.



Comfortable, Ultra-Breathable, Non-deformation
Multi-functional adjustment, 5° tilted forward, to maintain dentist a comfortable natural position when in surgery

Aluminium alloy base, stable and steady


360° switch silent casters, convenience and quiet.

Extraordinary Quality



Personalization

Optional
More options, More Functions



Built-in Scaler (Digital Display)



Built-in Micromotor (Digital Display)

Style
Customized



Standard



Top mounted



Trolley style

Quality assurance



Color
Breathable seamless microfiber leather
Choose your favorite



FS01
Modern brown



FS02
Coloured glaze blue



FS03
Gray



FS04
Begonia Red

Recommend color
FS02 Coloured glaze blue



FS05
Hermes orange



FS06
Olive-green



FS07
Tiffany Blue



FS08
light green